

EXPO in Rigaku software

Akihiro Himeda

Rigaku Corporation

10th, Sep, 2019

- Rigaku Powder XRD software overview
- Application example
 - Structure determination for a sample including impurities
- Structure refinement and validation
- SmartLab Studio II, Powder XRD plugin demo

- **Rigaku Powder XRD software overview**
- Application example
 - Structure determination for a sample including impurities
- Structure refinement and validation
- SmartLab Studio II, Powder XRD plugin demo

SmartLab Studio (-2016)

Bridge software which launches individual executables listed below.

- **PDXL**
 - Powder diffraction
- GlobalFit
 - XRR
 - High resolution rocking curve
- Nano-Solver (SAXS)
- 3D Explore
 - RSM & PF
- 2DP



On this system, multiple executables should be launched and closed by users.

SmartLab Studio II (2016-)

Integrated platform including measurement and analyses for X-ray.

Following functionality are provided as plugins.

- XRD Measurement
- **Powder XRD** (successor to PDXL)
- XRR (reflectivity)
- HRXRD
 - High resolution rocking curve & RSM
- Stress
- Texture
 - Pole figure & ODF
- SAXS


The screenshot displays the SmartLab Studio II software interface. The main window shows a 'Peak Profile View' with a plot of Intensity (cps) versus 2θ (degrees). The plot shows a series of peaks, with the most prominent one at approximately 18 degrees. The residual error is shown as a purple line at the bottom of the plot. The software is running a 'Structure determination' task for the dataset 'RulenicAcid_salicylicAcid'. The 'Phase Information' panel shows the phase name as 'Unknown' and the chemical formula as 'C₁₄H₁₀O₄'. The 'Candidate phase' panel shows a search result for 'RulenicAcid_salicylicAcid' with a chemical formula of 'C₁₄H₁₀O₄' and a concentration of 0.00 wt%.

SmartLab Studio II
Guidance™ expert system software
Version 4.4.0.0
Loading MachineLearningService v1.0 ...

Rigaku
© 2014 Rigaku Corporation. All rights reserved.

Powder XRD plugin – successor to PDXL

- Automatic peak processing
- Integrated qualitative analysis
 - Search/Match
 - Indexing
 - CIF
- Quantitative analysis methods
 - RIR quantification
 - Quantification using calibration curve
 - Quantification using Rietveld method
- Size & Strain
 - Scherrer's method and conventional Williamson-Hall plot
 - Fundamental parameter method: individual peaks and whole pattern
- Structure determination
 - Charge flipping, Direct space methods (parallel tempering method as in FOX)



Powder XRD plugin – Support for academic software

- EXPO
 - Direct methods
 - Direct space methods (SA)
 - Space group determination
- Import / Export Whole powder pattern fitting results
 - GSAS
 - FullProf
 - FOX
 - Rietan-FP
- Indexing
 - DICVOL06
 - ITO
 - N-TREOR
- Crystallography Open Database


Powder XRD plugin – Supported commercial databases

- ICDD
 - Phase identification (Search/Match)
 - Structure information for Rietveld refinement.
- Cambridge database
 - Restraint settings for Rietveld refinement


- Collaboration started on 2008 with EXPO2008
- **PDXL** Structure analysis package was released on Oct./2009 (PDXL ver. 1.5) together with EXPO2009.
- PDXL ver. 2.0 was released on 2011. Structure analysis guidance functionality was included.
- PDXL ver. 2.1 was released on 2013 together with EXPO2011.
- **Powder XRD** plugin Structure determination package was released on June/2018 together with EXPO2014.

- Rigaku Powder XRD software overview
- **Application example**
 - **Structure determination for a sample including impurities**
- Structure refinement and validation
- SmartLab Studio II, Powder XRD plugin demo


Co-crystal of flufenamic acid and salicylic acid




- Collaborative research with Chugai pharmaceutical company




- Measurement was performed on Rigaku SmartLab
- Sample was filled in a capillary and was rotated on the capillary spinner during measurement







- First attempt of indexing for cocrystal was failed
- Components (API* or co-former) were suspected to be included => tried to import their CIF files and found salicylic acid was contained
- Indexing was succeeded by removing impurity peaks.

*API: active pharmaceutical ingredient




- Rigaku Powder XRD software overview
- Application example
 - Structure determination for a sample including impurities
- **Structure refinement and validation**
- SmartLab Studio II, Powder XRD plugin demo


- Structure refinement with Powder diffraction needs to apply **restraints** to keep its structure reasonable.
 - > Reasonable partial structures are required.
- Mogul software from CCDC [Cambridge Crystallographic Data Centre] provides statistics for similar partial structures.
 - > It can be used for the target values and their restraint weight of bond distances and angles.



- Hydrogen bonds and Short contacts



- Void display empty space where one hydrogen can be included



- CIF export and launch checkCIF web page.

```

cif_core.dic 2.4.2 ftp://ftp.iucr.org/pub/cif_core.dic
cif_pd.dic 1.0.1 ftp://ftp.iucr.org/pub/cif_pd.dic
=====
# PROCESSING SUMMARY (IUCr Office Use Only)
=====
_journal_date_recd_electronic ?
_journal_date_from_coeditor ?
_journal_date_accepted ?
_journal_coeditor_code ?
=====
# SUBMISSION DETAILS
=====
_publ_contact_author_name 'ENTER NAME'
_publ_contact_author_address
;
ENTER ADDRESS
;
_publ_contact_author_fax 'ENTER FAX NUMBER'
_publ_contact_author_email 'ENTER EMAIL ADDRESS'
_publ_contact_author_phone 'ENTER PHONE NUMBER'
_publ_contact_letter
;
ENTER TEXT OF LETTER
;
_publ_requested_journal 'ENTER JOURNAL NAME HERE'
_publ_requested_category ? # CHOOSE FA FI FO FM EI EO EM AD or SC
_publ_requested_coeditor_name ?
=====
# TITLE AND AUTHOR LIST
=====
_publ_section_title
;
ENTER SECTION TITLE
;
;
_publ_section_title_footnote

```

checkCIF/PLATON (full publication check)

You have not supplied any structure factors. As a result the full set of tests cannot be run.

No syntax errors found.
Please wait while processing

CIF dictionary
Interpreting this report

Datablock: I

```

Bond precision:          = 0.0000 A
Cell:                   a=7.7875(10)   b=8.0610(11)   c=16.779(2)
                        alpha=93.2591(10) beta=91.0681(15) gamma=115.3082(11)
Temperature: 298 K
Wavelength=1.54059
=====

```

| | Calculated | Reported |
|--------------------|---------------------------|--------------------|
| Volume | 947.2(2) | 947.1(2) |
| Space group | P -1 | P -1 |
| Hall group | -P 1 | ? |
| Moiety formula | C14 H10 F3 N O2, C7 H6 O3 | ? |
| Sum formula | C21 H16 F3 N O5 | C42 F6 H32 N2 O10 |
| Mr | 419.35 | 838.70 |
| Dx, g cm-3 | 1.470 | 1.470 |
| Z | 2 | 1 |
| Mu (mm-1) | 1.074 | 0.000 |
| F000 | 432.0 | 432.0 |
| F000' | 433.70 | |
| h,k,lmax | 7, 8, 16 | 0, 0, 0 |
| Nref | 1950 | 1950 |
| Tmin, Tmax | | |
| Tmin' | | |
| Correction method= | Not given | |
| Data completeness= | 1.000 | Theta(max)= 50.000 |
| R(reflections)= | | wR2(reflections)= |
| S = | Npar= | |

The following ALERTS were generated. Each ALERT has the format **test-name_ALERT-alert-type_alert-level**. Click on the hyperlinks for more details of the test.

- Rigaku Powder XRD software overview
- Application example
 - Structure determination for a sample including impurities
- Structure refinement and validation
- **SmartLab Studio II, Powder XRD plugin demo**

Thank you for your attention!